





les encores®

THREADS

Product	Effect	Description	Recommended Treatment Areas
 <p>Mono</p>	<ul style="list-style-type: none"> Rejuvenation Tightening Smoothing 	<ul style="list-style-type: none"> Plain mono suture USP 5 Ultra-thin wall & sharp needle Allows for precise correction of fine lines Available in PDO, PLLA & PCL For beginner and advanced injectors 	<ul style="list-style-type: none"> Face & body skin firming Face shape & body contouring Any area requiring rejuvenation Any area of repair for fine lines & wrinkles
 <p>Screw</p>	<ul style="list-style-type: none"> Rejuvenation Tightening Smoothing Light volumization Reinforcement 	<ul style="list-style-type: none"> Coiled mono suture USP 5 ensures a higher concentration of material leading to higher collagen production Ultra-thin wall & sharp needle Allows for precise volumization of fine lines Available in PDO, PLLA & PCL For beginner and advanced injectors 	<ul style="list-style-type: none"> Face & body skin volumization Correction of linear wrinkles by volumizing Smoothing glabella Cellulite and scar repair Shallow marionette lines correction
 <p>HA</p>	<ul style="list-style-type: none"> Rejuvenation Tightening Smoothing Moderate volumization 	<ul style="list-style-type: none"> Plain mono or screw suture USP 5 Ultra-thin wall & sharp needle Allows for precise correction of fine lines Available in PCL For beginner and advanced injectors 	<ul style="list-style-type: none"> Face & body skin firming Face shape & body contouring Any area requiring rejuvenation and volumization Any area of repair for fine lines & wrinkles
 <p>Multi</p>	<ul style="list-style-type: none"> Rejuvenation Tightening Smoothing Heavy volumization 	<ul style="list-style-type: none"> 8-10 Mono USP 5 threads folded in half to create a fan of 16-20 threads inserted at one time. Ultra-thin wall & blunt tip cannula Available in PDO For beginner and advanced injectors 	<ul style="list-style-type: none"> Marionette lines and/or nasolabial folds Shallow Temples Body contouring and cellulite repair Jawline definition Chin enhancement
 <p>Rhino</p>	<ul style="list-style-type: none"> Volumization in bridge of the nose. 	<ul style="list-style-type: none"> Barbed suture USP 2 Ultra-thin wall & blunt tip cannula Available in PDO and PCL For advanced injectors 	<ul style="list-style-type: none"> Decreasing the width of the base of the nose Increasing the projection of the nose Adding volume to the bridge of the nose
 <p>Nasal</p>	<ul style="list-style-type: none"> Subtle lifting of the tip of the nose. 	<ul style="list-style-type: none"> Barbed suture USP 2 Ultra-thin wall & blunt tip cannula Available in PDO For advanced injectors 	<ul style="list-style-type: none"> Aligning the dorsum Aligning of the nasal wings Reducing the height of nasal wings
 <p>Eye/Lip</p>	<ul style="list-style-type: none"> Rejuvenation Tightening Smoothing Light volumization Reinforcement 	<ul style="list-style-type: none"> Mono & Screw suture USP 7 Ultra-thin wall & blunt tip cannula Designed for the delicate area around the eyes with minimal bruising. Available in PDO, PLLA and PCL For beginner and advanced injectors 	<ul style="list-style-type: none"> Tear trough correct Dark under eye circles Crows feet Correction of linear wrinkles Smoker's Lines Lip Enhancement

Product	Effect	Description	Recommended Treatment Areas
	<ul style="list-style-type: none"> Scaffolding Mild tissue repositioning Reinforcement 	<ul style="list-style-type: none"> Bidirectional, 4-dimensional oriented COGS USP 2 hold the repositioned skin in place Strong skin fixation with COGS in all directions Ultra-thin wall & ½ blunt tip cannula allowing for multiple placement attempts. Available in PDO & PCL For advanced Injectors 	<ul style="list-style-type: none"> Reinforcement and repositioning of the middle and lower third of the face Correction of gravitational ptosis of the eyebrows or neck
	<ul style="list-style-type: none"> Scaffolding Mild tissue repositioning Reinforcement 	<ul style="list-style-type: none"> Bidirectional, 4-dimensional oriented COGS USP 2 hold the repositioned skin in place Strong skin fixation with COGS in all directions Ultra-thin wall & blunt tip cannula Available in PDO & PCL For skilled and advanced Injectors 	<ul style="list-style-type: none"> Reinforcement and repositioning of the middle and lower third of the face Correction of gravitational ptosis of the eyebrows or neck
	<ul style="list-style-type: none"> Maximal fixation Dramatic tissue repositioning 	<ul style="list-style-type: none"> Ultra-thin wall & ½ blunt tip cannula allowing for multiple placement attempts. For use on patients with thick or heavy skin Available in PDO For advanced injectors 	<ul style="list-style-type: none"> Middle and lower face lifting & contouring Body lifting & contouring Correction of gravitational ptosis of the eyebrows or neck
	<ul style="list-style-type: none"> Ultimate fixation Maximum tissue repositioning 	<ul style="list-style-type: none"> Ultra-thin wall & ½ blunt tip cannula allowing for multiple placement attempts. For use on patients with thick or heavy skin Available in PDO For advanced injectors 	<ul style="list-style-type: none"> Middle and lower face lifting & contouring Body lifting & contouring Correction of gravitational ptosis of the eyebrows or neck

What makes Les Encres special?

Everything we do is to enhance your vision of beauty. The beauty of growing your business, helping your patients and setting the standard of excellence. We not only offer unmatched training material and business development tools specifically designed to solidify your return on investment, we also provide our partners with:

- A dedicated Clinical Specialist who is available to you to answer all your questions and ensure your team is confident and skilled on the art of thread insertion.
- A Custom Designed Thread Pouch. We designed our product pouch to resist heat and moisture, maintaining product stability and viability more so than any other thread pouch in the market today. This is a crucial step in quality product delivery. If the contents of the pouches are not diligently protected from heat, light and moisture the chemical compound of the sutures can begin to break down and oxidize.
- Peace of mind. Les Encres is registered with the FDA (#10058431) and offers a General Liability and Product Liability Insurance Policy. We are the only thread company in the industry to have \$1,000,000 liability policies, both General and Product liability.
- An Unsurpassed Thread Selection; PDO, PLLA PCL and HA PCL. Les Encres has THE largest thread selection in the market which allows you to fully customize your thread treatments and better meet your patient's needs.